


I d e a s t o S o l u t i o n s

**Elcam  
Medical**



# SWABBABLE STOPCOCK

Luer Activated Valve Integrated Stopcock



[www.elcam-medical.com](http://www.elcam-medical.com)

# SWABBABLE STOPCOCK

Luer Activated Valve Integrated Stopcock

A stopcock integrated with a luer activated Swabbable valve (or valves) allowing safe, needleless use of the stopcock ports. This offers leak-free aseptic fluid administration and sampling in the operating room, while keeping the IV line closed, preventing contamination of the I.V line.

The Closed Stopcock is cost effective as it eliminates the need for caps, "Hep-locks", cannulas or any other add-on components.



## Luer Activated Valve:

- Needle free connection - provides safety for health care personal and patients by preventing needle-stick injuries and contamination.
- Automatic closure after usage, preventing exposure of the fluid path to the atmosphere.
- Swabbable surface enabling aseptic use

## Applications

Suitable for all IV lines and anesthesia applications:

- Multilumen CVC
- "Closed system" infusion lines
- Blood sampling and infusion

## Technical specifications

### • Standard Luer Connections

Easy assembly and integration into kits and set, suitable for luer-lock and luer-slip male connectors.

### • Minimal Residual Volume

Residual volume at valve port is 0.12cc

### • High Flow Rate

Minimum flow rate of 150 ml/min under 1.5 psi.

### • Other

Lipid Resistant

Latex Free

Non-hemolytic

## Available Configurations

### • High Flow (large bore)

- One Port
- Two Ports

### • Standard Flow

- One Port



Head Office / Manufacturer  
**Elcam Medical**  
BarAm 13860  
Israel  
Tel: 972-4-698-8120/1/2  
Fax: 972-4-698-0777  
sales@elcam.co.il  
www.elcam-medical.com

USA  
**Elcam Medical INC.**  
Suite 620, 2 University Plaza,  
07601 Hackensack, NJ  
Tel: +1 201 4571120  
Fax: +1 201 4571125  
info@elcam-medical.com  
www.elcam-medical.com

Europe  
**MedNet GmbH**  
Berkstrasse 10  
48163 Muenster Germany  
Tel: +49 251 32266-0  
Fax: +49 251 32266-22  
info@medneteuropa.com  
www.medneteuropa.com

China  
**Sheng Enterprises Ltd.**  
Suite 2003 Universal Mansion  
No.168 Yu Yuan Rd.  
SHANGHAI, 200040  
Tel: 86-21-62497255/6/7  
Fax: 86-21-62497259  
johnliu@sheng-cn.com



ISO 13485  
ISO 9001